THE RADIX LENS

A CLEAR CHOICE IN PASSIVE TECHNOLOGY

As the popularity of computer-assisted therapy continues to grow, applications continue to be introduced that place new demands on the technology. Optical tracking is now being used in environments that require markers to be robust to liquid and tissue contamination. They also need to be easily cleaned so they can quickly recover tracking.

The NDI Radix™ Lens has been designed to respond to these new demands. The lens has a smooth plastic surface that naturally sheds liquids and is easy to wipe to recover tracking. The Radix Lens also has the unique characteristic of being inherently resistant to partial occlusions, further enhancing its robust performance characteristics.

Designed for use with the Polaris™ and Polaris™ Spectra systems, the Radix Lens may be the right choice for your next application.

Features

- Facilitates accurate optical tracking in harsh environments
- Solid contamination is easy to identify and remove resulting in more efficient tracking recovery
- Partial occlusions have minimal impact on tracking performance
- Minimizes interruption to surgical workflow
- Low profile marker eliminates marker-to-marker interference on planar, multi-marker tools

Partial Occlusion

The Radix Lens* has the unique characteristic of being inherently resistant to partial occlusions. For passive spheres partial occlusions generally cause a tool to stop tracking.

A partial occlusion covering one side of a Radix Lens will be reflected on the opposite side. Due to this symmetrical occlusion, there is little effect on the marker centroid and therefore tracking loss is less likely.

Tool Integration

In contrast to the NDI passive sphere, the Radix Lens is not used directly in the form that it is supplied. The Radix Lens must first be incorporated into a tool or mounting base so that it may be attached to a passive tool. Please contact NDI for details regarding integration of the Radix Lens with passive tools.


Figure - Mechanical Details

Off-Angle Viewing

The Radix lens has a maximum marker angle of $\pm 70^\circ$.

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