

# Accurate and Reliable

## World leading eye control unit



## Tobii CEye is an eye control unit that simply works



### **What is eye control?**

Eye control is an alternative and complementary input method to using a mouse and keyboard. By looking at a screen you control 'the mouse' and you click by blinking, dwelling (staring at the screen for a certain length of time) or using a switch.

When talking about eye tracking all components have to be considered – hardware, software, peripherals, service and support – and how well they work together. Designed and built by Tobii, the Tobii CEye eye tracker works flawlessly with compatible Tobii communication devices.

### **Easy to use**

It takes only minutes to set up the Tobii CEye. Users do not need to do or wear anything; they simply place themselves in front of the CEye-equipped device. After a fast, one-time calibration, no additional adjustments are needed.

## **Eye control that is simple to set up and use, accommodates large head movements, and is reliable and extremely accurate.**

Tobii eye control is fully automatic and so easy to use that people with all levels of cognitive development can use it.

### **Inclusive and reliable technology**

CEye can be relied upon to perform. It operates with a high level of precision in nearly all light conditions. It also works with most users, regardless of ethnic origin – and even glasses or contact lens users.

### **Freedom of movement**

The Tobii CEye has one of the market's largest track boxes, which refers to the size of the imaginary box within which you can move your head.

Thanks to the large track box it isn't critical that you remain in the same place. As long as you keep within the track box you can move freely without having to worry about interrupted tracking or diminished functionality.

And should you move outside the track box, it isn't a problem. The Tobii CEye will start tracking you again as soon as you reenter, without requiring the system to recalibrate.

This unique ability to deal with large head movements has opened up eye control for people who previously were unable to use other eye controlled systems, including children and users with Cerebral Palsy.

### **Powerful direct control**

Tobii devices equipped with CEye use a technique that allows you to interact directly

with objects on the screen. This approach decreases stress and increases speed dramatically. It allows you to be much more efficient for longer periods of time than ever before.

### **Communication software of choice**

Leading industry software has been Tobii-enabled and modified to be more effectively accessible through eye control. As a result you can choose the best communication package for your needs, so you get the most out of your Tobii device.

### **Third party program**

With the Tobii third party program developers can more easily design applications that work with the ever increasing number of Tobii units.

The program ensures a steady flow of software applications that help Tobii users to communicate, learn, work and play using eye and direct control.

© Tobii®. Illustrations and specifications do not necessarily apply to products and services offered in each local market. Technical specifications are subject to change without prior notice. All other trademarks are the property of their respective owners.

#### **NORTH AMERICA**

Tobii ATI  
333 Elm Street  
Dedham, MA 02026  
Toll-Free: 800-793-9227  
Phone: 781-461-8200  
Fax: 781-461-8213  
Web: [www.tobiiati.com](http://www.tobiiati.com)  
[sales.us@tobii.com](mailto:sales.us@tobii.com)

#### **EMEA**

Tobii Technology AB  
Karlsrovägen 2D  
S-182 53 Danderyd  
Sweden  
Phone: +46 8 663 69 90  
Fax: +46 8 30 14 00  
[sales@tobii.com](mailto:sales@tobii.com)

#### **ASIA**

Tobii Technology, Ltd.  
3-4-13 Takanawa, Minato-ku  
Tokyo 108-0074  
Japan  
Phone: +81-3-5793-3316  
Fax: +81-3-5793-3317  
[sales.jp@tobii.com](mailto:sales.jp@tobii.com)

#### **CENTRAL EUROPE**

Tobii Technology GmbH  
Niedenau 45  
D-60325 Frankfurt am Main  
Germany  
Phone: +49 69 24 75 03 40  
Fax: +49 69 24 75 03 429  
[sales.de@tobii.com](mailto:sales.de@tobii.com)